




Breathe Health, Breathe  
Molecular Hydrogen



[www.HydrogenTherapyH2E.com](http://www.HydrogenTherapyH2E.com)



## USER MANUAL

### Hydro Gen Version 2 1200ml

#### Hydrogen oxygen machine

(Electrolytic water hydrogen and oxygen production equipment)



# Contents

1. Notes for Use (Please read)	3-5
2 Technical parameters	6-6
3 Environmental conditions of the operation place	6-6
4 Instructions for Use	7-8
5 Product introduction and operation instructions	8-14
6 Maintenance and cleaning	15-15
7 Transport and Storage	16-16
8 Common faults	16-18
9 Pollution Control Labels of Electronic Information Products	18-19
10 Liability of Manufacturer	20-20
11 Packing List	21-21
Warranty card	

● Please read  
● Warm pro

No. 7




1

Notes for Use (Please read)


In order to ensure the safe and effective use of Hydrogen Health Care Machine in good condition, please read the following notes carefully first, and then read the operating instructions of the product carefully.

⚠


WARNING




When using, there shall be no ignition device, or action generating static electricity to prevent explosion!




The collection of hydrogen and oxygen is prohibited; Other uses are prohibited.




When this product is in use, smoking is prohibited, and it shall be kept away from open fire!




No sundries and oil containers shall be placed on the top of the product.




When using, it is forbidden to tilt, invert and recline this product.



Please use distilled water or purified water, and please add proper amount distilled water or purified water before starting it for use!




When this product is used by children or persons with impaired behavior, there shall be guardians to attend on them;




Please use distilled or purified water. Please add an appropriate amount of distilled or purified water before turning it on.

⚠


ILLUSTRATION




There is something to remind you「caution」「danger」「warning」



Compulsory or specified sexual acts




Prohibited acts




Disaggregation prohibited




⊘ Do not use in the following environments




High temperature and humidity environment




Oily environment




Bathroom and other water environment




Water vapor rich environment



Dusty environment




Flammable and explosive (gas or dust) environment




Environment with strong electromagnetic interference


⊘ Please observe the following prohibitions, otherwise it may cause failure or danger.




Do not splash water or other liquid onto the body




Do not force, twist or tie the power cord




Do not touch and plug with wet hands




In case of lightning, power on shall be terminated and power plug shall be unplugged




Do not use sharp objects to touch and scratch the screen



Do not share a socket with other equipment



Do not place heavy objects on the power line



No heavy blows

⚠

- Avoid generating any sparks near this product during use, including sparks caused by various frictional static electricity;
- Children over eight years old and persons with reduced physical, sensory or mental abilities or lack of experience and knowledge should use this product under the supervision or guidance of a qualified person.



- Minors and children are not allowed to operate, play, or use this product;
- Product cleaning and maintenance must not be performed by children without supervision;
- The use and purpose of the machine cannot be changed, and hydrogen and oxygen cannot be collected for breathing or other uses (for example, the use of oxygen bags is prohibited).

- Please use it in a ventilated environment;
- Power shall be cut off after each use (turn off the main power supply);
- When this product is not in use, please turn off the main switch and unplug the power plug.
- When moving the main engine, please do not place it in the place with abnormal high temperature, or it may cause breakdown and other hazards.
- This product shall not be placed in the following environment: near the heat source and open and dying fire, wet, unobstructed, with smoke and pollution, too high or too low temperature;
- Please try to avoid the places with bad effects such as moisture, dust, salt, sulfur, etc., or it may cause fire, breakdown and other hazards.



- This appliance must not be used by persons (including children) with reduced physical, sensory or mental capabilities or lack of experience and knowledge unless they have been given supervision or instruction;
- Please keep the nasal oxygen tube out of the reach of children to prevent the connecting tube from getting entangled in the neck and causing danger;
- When this product is used by individuals at home, it should be used under the guidance of professionals;
- Do not use it with other therapeutic equipment at the same time;
- Do not place ME EQUIPMENT in a location where it is difficult to operate the disconnect device, as this may result in a safety risk; please read the instructions for use carefully and obey all instructions and warnings.



- Please do not disassemble the machine. Non-professionals or authorized personnel of the company are not allowed to disassemble, modify or repair this product.

- Regularly wipe off the dust on the host and power plug. If there is dust on the plug, moisture may cause poor insulation and lead to fire;
- Do not use it in an environment with electromagnetic interference;
- The wires and power supply head must not be damaged. Use after damage may cause electric shock, short circuit, or fire. For repairs of the wires and power supply head, please contact the service hotline.

2

Technical parameters

Hydrogen oxygen machine (Electrolytic water hydrogen and oxygen production equipment)				
Specification and model	YSR-450S	YSR-900S	YSR-1200S	YSR-1500S
Mixture ratio of hydrogen and oxygen	Hydrogen 66.6%, Oxygen 33.3%			
Output flow	300/450 mL/min	450/900 mL/min	450/900/1200 mL/min	450/900/1500 mL/min
Rated voltage	100-240V~ 50/60Hz			
Rated input power	160W	260W	300W	360W
Electrolytic cell technology	Hydrogen Production by Electrolytic Water with SPE Proton Exchange Membrane			
Hydrogen production water	Distilled water (purified water)			
External Dimension (L*W*H)	218x422x419 (mm)			
Complete Machine Net Weight	9.0 (kg)	9.2 (kg)	9.5 (kg)	9.8 (kg)

3

Environmental conditions of the operation place

1. Environment temperature: 5°C~40°C;
2. Relative humidity: ≤ 90%;
3. Atmospheric pressure: 86 kPa~106 kPa;
4. Power supply voltage: AC100-240V, 50/60Hz;
5. The use of distilled water or purified water shall meet the requirements of GB17323-1998;
6. No direct sunlight, and no open fire;
7. No large dust and conductive ions, and no acid, alkali and other corrosive gases;



- Note: When the storage temperature is below 5°C, the product shall be placed in the environment with normal working temperature for more than 4 hours before use.
- When the environment temperature is 5 ~ 10 °C, the machine will be in the preheating mode when it is first started, and the output flow rate of the preheating mode is low.



## 4 Instructions for Use

 Please note the following.

1. Before starting operation for the first time, please open the filler cover to add distilled water or purified water and make sure to meet the water level requirement of the water tank in time. Wait for 15 min after adding the water before starting operation.

15min

2. During use, please use the water tank observation window on the right to confirm that the water volume in the water tank should be controlled between the minimum water level (Min) 1600mL and the maximum water level (Max) 2500mL. Do not exceed the maximum water level;  
3. During use, if the water level in the water tank is lower than the lowest water level, the water shortage indicator light on the host will light up/voice will prompt. At this time, water needs to be added to the normal water level range to ensure the normal operation of the product;

water  
level  
(MIN)  
(MAX)

4. Water quality requirements: It is required to use pure water or distilled water for electrolysis. When the water quality TDS is greater than the set value, an alarm is triggered; hard ions (such as calcium and magnesium ions) in unqualified water will deposit, block the electrolytic cell, and cause damage to the electrolytic cell;  
Note: If the machine uses non-specified water, it will not be able to enjoy the three-guarantee service.

Water  
quality

5. When the product is out of service for a long time or needs to be transported, it is necessary to discharge the water from the water tank.  
6. During the process of using the product, if the water source that does not meet the requirement is used, it can be used normally only by carrying out the function of rinse drainage many times, adding new water that meets the requirement, and ensuring the water quality monitoring is normal.  
7. When the product is in use, the electrolytic cell and water tank shall be kept clean, since the long-term use of the product may lead to the breed of microorganism, resulting in turbid water and affecting the hydrogen oxygen production and service life of the product. Rinse Drainage shall be carried out according to the prompt of water quality monitoring function during the use of the product.  
8. It is recommended to change the water in the water tank every 1 to 2 weeks (use new water).

Cleaning  
and  
drainage  
function

9. The small holes on the filler cover of the water tank of the product are the outlets of the water tank, which cannot be blocked, and the filler cover of the water tank cannot be arbitrarily replaced!  
10. During the process of transportation, water shall not be stored in the water tank to prevent water spillage from damaging the electrical components;  
11. The filter element has a lifespan of one year and needs to be purchased separately. The filter element purifies the water quality to better protect the water used by the machine;

Water  
tank

12. During the use of the product, the two sides and the rear side of the main engine shall be kept above 15cm with other objects, and the cooling window shall not be covered, so as to prevent the blockage causing too high temperature and causing downtime or product failure;

Kept  
above  
15cm

13. If any adverse reactions occur or appear during the use of this product, please immediately discontinue the use and contact the supplier of the product or doctor.


Contact

## 5 Product Introduction and Operation Instructions

### 1. Unpacking


- 1) Please handle with care when unpacking and transporting;
- 2) Check the carton or other packaging for obvious damage. If there is any damage, please notify the transporter or dealer in time;
- 3) Take out the product and accessories and check whether the accessories are complete according to the packing list;
- 4) Check the outer surface of this product for damage such as chips and dents. If there is any damage or incomplete accessories, please notify the transporter or dealer in time.

### 2. Diagram of the complete machine







Back of host



Top of host



Condenser connection Schematic diagram (optional)





Product Features: Reduce and prevent water vapor from condensing into water droplets in the nasal oxygen tube;  
Instructions: For details, please refer to the usage instructions on the condenser packaging;  
Remark: This condenser is an optional, non-standard accessory. If necessary, please contact the dealer to purchase it.

3. Start-up instructions

1. Connect the power supply


- 1) Take out the power cord, connect the power supply, and press the power switch to "I". The power indicator light will light up and a voice prompt will sound at the same time.
- 2) The display screen lights up: the time set when it was last started is regularly displayed, and the accumulated running time of the machine is displayed at the same time; the current water quality TDS value is displayed, the default voice is on, and the running icon is always on.



2. Fill the water tank with water

- 1) Unscrew the water tank cover on the top of the host and add distilled or purified water; confirm that the water level requirements of the water tank are met through the water tank observation window on the right side of the product. The water tank is marked with an upper and lower water level limit, which needs to be controlled from the lower water level limit (Min) to the water level. Between the upper limit (Max), please do not exceed the maximum water level!

Water level diagram








2. Fill the water tank with water



- When the user fills the water tank with water, provided that the device is powered on and the water level does not reach the minimum level, the host will alarm and the system will give a voice reminder: "Water shortage reminder, please check the water level of the water tank and add an appropriate amount of distilled or purified water.", please add it between the lower limit of the water level (min) and the upper water level (max) of the water tank, and do not exceed the maximum water level;

Note: The silicone plug is only used to prevent water overflow during transportation and is not used normally.

- First remove the silicone plug, and then place the purification filter element (the filter element should be replaced once a year);
- If the water filling in the water tank exceeds the maximum water level, drainage needs to be performed; pull out the drainage pipe at the rear of the main unit, open the silicone drainage plug, and drain the excess water until the required water level is reached. Close the drainage plug to avoid water overflow and causing the main unit to malfunction. abnormal;
- The water tank used for machine operation is recommended to be replaced every 1 to 2 weeks.



Hold the handle of the hydrogen production observation window with one hand, and steady the machine with the other hand to remove the hydrogen production observation window (see Figure 1 below); open the top plastic cover and pour in an appropriate amount of distilled or purified water to ensure that the water level does not exceed the marked maximum water level. ;Tightly cover the plastic cover and ensure it is sealed; then connect the observation window to the main unit and confirm that it is installed in place to prevent loosening and air leakage.




Note: It is recommended to change the water in the observation window once a week.

3. Hydrogen production observation window

4. First boot




Note: For first time use, add pure water to the upper middle of the water level in the water tank, and test the drainage after soaking for half an hour;

Test drainage: Open the drainage hole pipe at the rear of the host, refer to the "Water Tank Filling" operating instructions, and drain a small amount of water and air to ensure that the water path of the electrolytic tank inside the host is smooth.



5. Connect the water collector and nasal oxygen tube

- Take out the nasal oxygen tube and water collector, install the water collector on the water collector holder, then connect one end of the water collector connecting pipe to the air inlet (thick) on the water collector, and the other end to On the hydrogen and oxygen gas outlet on the main unit (connect according to the arrow instructions);
- Connect the connecting end of the nasal oxygen tube to the air outlet (thinner) on the water collector. The nasal oxygen tube is worn by the user to inhale hydrogen;



**It is prohibited for multiple people to use it at the same time!**

- Please connect the water collector before each use of the product. Check and clean the water accumulated in the water collector or water condensed from water vapor during use of the product. This can prevent the nasal oxygen tube from being inhaled into the nasal cavity during breathing and causing discomfort;
- During the use of the machine, the nasal oxygen tube used for the hydrogen and oxygen outlet must not be folded or blocked, which may easily lead to poor air outlet and damage to the machine.



- The nasal oxygen cannula is a consumable accessory. For personal hygiene, please clean and disinfect it before use;
- It is recommended to replace the nasal oxygen tube regularly every two months. If it is used frequently, it is recommended to replace it every month;
- It is normal for the nasal oxygen tube to accumulate water vapor during use. You can shake the nasal oxygen tube before use to remove the water droplets or water vapor.



## 6. Boot and run

After pressing the hydrogen-oxygen key to start, check whether there are bubbles in the hydrogen production observation window. The presence of bubbles indicates that the product is working normally. Users can use the product correctly according to the operating instructions in the manual, or according to the knowledge and recommendations of professionals.



- 1) For the first time water filling and startup or when it has not been used for a long time, you need to wait 15 to 30 minutes after filling the water tank before running the hydrogen production function to ensure that the electrolyzer is not short of water;
- 2) The water tank should be kept clean for a long time. Please clean, drain and replace the designated water according to the water quality requirements.

## 4. Button operation instructions

★ After unpacking, you need to wait 15 to 30 minutes after the first water injection, and then start operation again to ensure that the electrolyzer is not short of water.



- 1 Press this button, YSR-450S model sets the hydrogen and oxygen flow rate to switch between 300/450; YSR-900S model sets the hydrogen and oxygen flow rate to switch between 450/900; YSR-1200S model sets the hydrogen and oxygen flow rate to switch between 450/900/1200; YSR-1500S model sets the hydrogen and oxygen flow to switch between 450/900/1500; at the same time, the LED screen It should be displayed that when the voice is turned on, every time this button is pressed, the voice will broadcast the corresponding set traffic;
- 2 Press and hold this key to switch between "night mode" and "normal mode"; "night mode": screen indicator The brightness of the indicator light is dimmed, and the ambient lights of the humidification cup and water tank are turned off.



Press this button to set the selectable time, which can be switched cyclically between 30 minutes, 40 minutes, 50 minutes and 1~8 hours in steps of 1h. It is also equipped with a voice prompt function; a time memory function that records every time it is turned on. The time set when the run was last started.



Press this button to switch between voice on and voice off. The corresponding voice icon on the LED screen will light up and go out. When lit, the voice is on and has a voice prompt function.



Press this button to start and stop the running function, and there will be a voice prompt at the same time. When the hydrogen and oxygen mixture is turned on, the running icon starts to flash, and the "Hydrogen 66.6% Oxygen 33.3%" icon lights up at the same time; the hydrogen and oxygen mixture function stops. , the running icon is always on, and the "Hydrogen 66.6% Oxygen 33.3%" icon is extinguished at the same time.


## Grand total time

The output unit is 1 hour, which means the current displayed time is the total accumulated time in the running state of the machine.

## Water quality show

- 1 In order to ensure the water injection requirements of the machine, the machine has a built-in water quality monitoring sensor;
- 2 Displays the TDS value of the current water quality of the machine. Values from 0 to 14 are normal. The higher the TDS value, the more dissolved substances are contained in the water;
- 3 Usually refers to the concentration of Ca<sup>2+</sup>, Mg<sup>2+</sup>, Na<sup>+</sup>, K<sup>+</sup> and other ions in the water, which means the water is impure;
- 4 When the water added by the machine does not meet the requirements or after a period of use, the water quality changes, the machine will trigger the water quality alarm E03 to indicate that the water quality is not up to standard, and require the unclean water to be discharged and pure water to be added again.

## 5. Power socket description



- 1) Drainage hole: When draining the host water tank, you need to pull out the drain pipe and remove the drain plug to drain. Pay attention to the containment and treatment of the drainage. Do not drink;
- 2) Power switch, O closes, I opens;
- 3) Fuse: fuse holder, equipped with F5AL250V (5\*20mm) fuse tube;
- 4) Power socket: Connect the product to AC100-240V power supply.



How to replace the fuse:

- 1) Turn off the power and unplug the power cord from the power socket;
- 2) Use a screwdriver to pry open the fuse tube holder under the power socket and take out the fuse tube;
- 3) Replace the fuse specification to F5AL250V; if the fuse burns out again after replacing it, please contact our company's maintenance center.



6

Care, cleaning and maintenance

1. Check the status of the product during use

- 1) When using the product, you should pay attention to the following status checks:  
\* The displayed content should be correct;  
\* Key operations should be normal;  
\* The water levels in the water tank and observation window should be normal during use.
- 2) If any abnormalities are found in the above checks, the following handling methods can be adopted:  
\* Stop using it and contact your local dealer;  
\*Consult our company’ s after-sales service center.

2. Daily maintenance and cleaning

- 1) The power cord should not be deteriorated, deformed or damaged;
- 2) The water in the water tank and the observation window should be cleaned and replaced frequently to ensure that the product meets the water quality requirements; the water in the water tank is recommended to be cleaned and replaced once every 1 to 2 weeks, and the water in the observation window is recommended to be cleaned and replaced once a week;
- 3) When cleaning the observation window, first rinse it with light detergent, then rinse it with clean water to ensure hydrogen hygiene;
- 4) When cleaning the observation window, do not touch the internal filter element to prevent it from falling off;
- 5) When it is not used for a long time, it is necessary to clean and drain the water in the water tank and observation window to prevent bacteria from growing if it is not used for a long time;
- 6) The nasal oxygen tube is a consumable accessory. For personal hygiene and long-term use, please disinfect it frequently or replace it regularly;
- 7) The host does not require professional maintenance. It only requires routine cleaning. Use a soft cloth to wipe off the dust on the surface of the host. Keep the host dry and do not allow liquid to penetrate into the host.
- 8) The heat dissipation windows and back panel on both sides of the host can be lightly brushed with a small brush to avoid accumulation of dust, which may affect the heat dissipation and use of the machine.

3. After not using it for a long time, check it before using it again

- 1) Please wipe off the dust or dirt on the cord, power plug, and main unit, and keep them clean;
- 2) Please confirm whether the buttons and indicator lights are functioning properly;
- 3) After completing the above checks, please follow the initial power-on steps. For details, please refer to the user manual.

7

Transportation and storage

Packaged products should be transported and stored under the following environmental conditions

- 1) Ambient temperature: 0°C~55°C;
- 2) Relative humidity: no more than 93%;
- 3) No corrosive gas and well-ventilated environment.



- 1) This product should be stored in a room without strong sunlight, corrosive gases and well-ventilated. During transportation, avoid severe shocks and inversion;
- 2) When the storage and transportation temperature is lower than 0°C, the product should be placed in a room temperature environment for 48 hours before use, and then operated according to the instructions to avoid being affected by the environmental temperature of storage and transportation and causing product failure.

8


Common malfunctions

When the equipment detects a fault, the machine will stop running, the LED screen will display the corresponding fault code at the scheduled time, the corresponding fault icon on the LED screen will flash, and there will also be a buzzer and voice broadcast prompt function.

Serial No.	Fault Phenomenon	Inspection the Reasons	Maintenance Guide
1	When the power supply is turned on, the buzzer of the mainframe has no sound, the screen is not bright, and the product can not run.	1-Poor contact between power cord plug and socket	Recheck the plug-in socket
		2-No power output from the socket	Move to a socket with power output
		3-Fuse blown	Replace the fuse
2	Hydrogen-oxygen outlet has no flow output	1-Air system leakage	If such failure is found, stop using immediately. It must be sent to the manufacturer for repair and replacement. Do not disassemble it by oneself! Otherwise, at your own risk!
		2-Poor sealing of air and water separator	
		3-The system power supply is damaged	

Serial No.	Fault Phenomenon	Inspection the Reasons	Maintenance Guide
3	Mainframe does not work, and alarm sounds; Display screen displays "E01"	It is triggered when the water level is below the minimum water level limit;	Water shortage reminder, please check the water level of the water tank and add proper amount of distilled water or purified water;
4	Mainframe does not work, alarm sound ; Display screen display "E02" (Not applicable to this model)	Air plug: it is triggered when the hydrogen pressure is greater than the set maximum pressure set;	Air plug remind, please check whether or not the accessories are folded or pressed, or blocked!
5	Mainframe does not work, alarm sound ; Display screen display "E03"	Water quality: it is triggered when water quality is not up to standard;	Water pollution reminder: the current water quality does not meet the requirements, please timely replace "distilled water" or purified water;
6	Mainframe does not work, alarm sound ; Display screen display "E04"	Overheat: it is triggered when the temperature exceeds the set value;	Check whether the cooling windows on both sides of the mainframe are blocked, and start it again after cooling;
7	Mainframe does not work, alarm sound ; Display screen display "E06"	Improper external operations lead to the product triggering the safety protection device, and please operate and use the product correctly.	Please follow the instructions and operate in correct way, and read precautions.
8	Mainframe does not work, alarm sound ; Display screen display "E07"		
9	Mainframe does not work, alarm sound ; Display screen display "E08"	The low temperature alarm is triggered due to low ambient temperature	Place of use Ambient temperature: 5°C~40°C
10	During the use of the product, there are many fogs or water droplets in the hydrogen-oxygen tube.	1-Due to the large difference between the internal temperature and the ambient temperature during the long-term operation of the mainframe, water mist and water droplets are produced; 2-The water collection filter is too full; 3-The water in the hydrogen production observation window is too full, which has exceeded the maximum water level line; or the water temperature of observation window is too high;	Remove the water collector filter, open the filter, and pour out the water stored inside. Discharge water droplet in the hydrogen-oxygen tube or dry it. Wait for the mainframe to rest for a period of time, and restart it;  Remove the hydrogen production observation window and pour out the excess water.

Serial No.	Fault Phenomenon	Inspection the Reasons	Maintenance Guide
11	If water is found to flow out of the gap at the bottom of the mainframe.	Silicone hose is aging, and waterway interface is loosening	Send it to the manufacturer for inspection, repair or replacement



If other faults cannot be eliminated, please contact the dealer or our company's service hotline.

Other faults:

- 1.Non-professionals or authorized personnel of the company are not allowed to disassemble, modify, or repair. If there is any other failure with this product, you need to contact the local dealer or the company's technical service department;
- 2.After tearing off the warranty label on the back of the console, free warranty will no longer be provided.

9

Electronic information product pollution control label

Component and Part name	Toxic and hazardous substance or element					
	Lead (Pb)	Mercury (Hg)	Cadmium (Cd)	Hexavalent chromium (Cr6+)	Polybrominated biphenyls (PBBs)	Poly brominated diphenyl ethers (PBDEs)
Plastic shell	0	0	0	0	0	0
Circuit board	X	0	0	0	0	0
Connecting line	0	0	0	0	0	0
Metal components	0	0	0	0	0	0
Label class	0	0	0	0	0	0



Note: circuit board includes printed circuit board and its components and parts, such as solder, resistors, capacitors, integrated circuits and connectors, etc. The connecting wire includes copper wire and its raw material composed of it, such as solder and insulating material, etc.

0 : it indicates that the content of the toxic and hazardous substance in all homogeneous materials of the components and parts is below the limit specified in the SJ/T11363-2006 standard.

X: It indicates that the content of the toxic and hazardous substance in at least one homogeneous material of the components and parts exceeds the limit requirements specified in SJ/T11363-2006 standard.

Note: for components and parts marked with "X" above, some contain hazardous substances beyond the standard due to the limitation of the current industry technical level, so it is unable to achieve reduction or replacement temporarily.

## 2. Environmental protection use period



The number in this logo indicates that the environmentally friendly use period of the host under normal use is 10 years.

## 3. Product expiry date instructions

This product is a health electronic product, and its service life is usually related to the following factors:

- 1) Usage environment, such as network power supply stability, temperature, humidity, dust, etc.;
- 2) Intensity of use, such as duration of each use, number of uses, intensity level of use, etc.;
- 3) Maintenance conditions, such as regular dust removal, etc.;
- 4) When used under normal conditions, the product has a service life of 10 years. If it still needs to be used after 10 years, it should be tested.



## 4. Environmental protection

- 1) Do not dispose of this device and its accessories as ordinary household waste;
- 2) Please comply with local regulations for the disposal of this device and its accessories and support recycling initiatives.

# 10

## Manufacturer's Responsibilities

- 1 If the products produced do not meet the standard requirements after identification, our
- 2 company will be responsible for repair, replacement and return.  
If any other property losses or personal injuries are caused by defects in the product, our company shall bear corresponding responsibilities. However, our company shall not be liable in any of the following circumstances:
  - 1) When the product is put into use, the defect causing damage does not yet exist.
  - 2) The purchaser and operator did not read the instruction manual carefully and did not use it in accordance with the instructions.
  - 3) Failure to comply with product installation, use, maintenance, transportation, storage, etc. regulations, resulting in product, personal, and property losses.
  - 4) Matters specified in the warranty card.
- 3 Manufacturing quality problems that occur within the specified quality warranty period will be repaired for free for the user and corresponding technical services will be provided. Please refer to the warranty card for the warranty period agreement.



11

Packing List

1. Unpacking precautions

- 1) When unpacking, the front of the packaging box should be facing up and not turned upside down;
- 2) Check the integrity of the outer packing box;
- 3) When unpacking and handling, please handle it with care without damaging the machine;
- 4) The user accepts the product according to the packing list.

2. Packing configuration list

Part Name	Quantity
Host	1 set
power cable	1 root
Nasal oxygen tube	2 root
water collector	1 set
Instructions (including warranty card)	1 serving
Certificate	1 serving
fuse	2 serving
Purification filter element	1 set

3. Optional accessories

Part Name	Quantity
Condenser	1 set

Note: 1.This list is the company’s factory standard configuration (excluding optional accessories) and is subject to change without prior notice.  
2. Water collectors and humidification cups are consumables and need to be replaced regularly; starting from the period of use, it is recommended to replace them at least once a year.